



CARR Australia Pty Ltd
ABI Distribution Pty Ltd

1 Spireton Place, Pendle Hill.
 PO Box 335, Pendle Hill.
 NSW 2145, Australia.
 Telephone 61-2-9636 9525
 Facsimile 61-2-9636 8844
 www.carrgroup.com.au

**APPROVED CLEANING
 PRODUCTS FOR MATTRESS
 COVERS, UPHOLSTERY
 TOPS AND SHOWER
 TROLLEY TUB FABRICS:**

CLEANING PRODUCT:	MAXIMUM CONCENTRATION:
Clinell Hand & Surface Sanitizer Disinfectant Clinell Universal Sanitising Wipe Clinell Detergent Wipes Clinell Sporidical Wipes (GlobalMedics, Australia & New Zealand)	Use according to the manufacturers specifications.
PAL Detergent Wipes (GlobalMedics, Australia & New Zealand)	Use according to the manufacturers specifications.
Mediwipes Disinfection (Sulco Limited, New Zealand)	Use according to the manufacturers specifications.
PAL Disinfectant Wipes (GlobalMedics, Australia & New Zealand)	Use according to the manufacturers specifications.
Surface Sanitising Wipes Neutral Detergent Wipes (Reynard Health Supplies, Australia, New Zealand & Europe)	Use according to the manufacturers specifications.
PureGreen 24 (Pure Green 24, America)	Use according to the manufacturers specifications.
Spring Clean (Eco Lab, Australia)	5% Spring Clean (1 part Spring Clean to 20 parts water 50ml per litre)
Ecopine (Worldwide)	5% Ecopine (1 part Ecopine to 20 parts water 50ml per litre)
Actimax (Worldwide)	Use according to the manufacturers specifications.
Viraclean (Whiteley Medical, Australia)	Use according to the manufacturers specifications.
Virkon (Dupont, Worldwide)	5% Virkon (1 part Virkon to 20 parts water 50ml per litre)
Trigene (Bio Strategy, Australia & New Zealand)	Use according to the manufacturers specifications.